VIDAS® FPSA (FPSA)

IVD

VIDAS® FPSA is an automated quantitative test for use on the VIDAS family instruments, for the quantitative measurement of the free fraction of prostate specific antigen (PSA) in human serum or plasma (lithium heparin or EDTA) using the ELFA technique (Enzyme Linked Fluorescent Assay).

SUMMARY AND EXPLANATION

Prostate-specific antigen (PSA) is a glycoprotein which belongs to the kallikrein family. PSA has a molecular weight of 30,000 Daltons (1). PSA is principally produced by the glandular epithelium of the prostate, and is secreted in the seminal fluid. PSA is also present in urine and blood. PSA acts on seminal fluid to fluidify and increase sperm mobility.

PSA is present in blood in three main forms (2, 3). The most important immunoreactive form is PSA bound to Alpha-1-antichymotrypsin (PSA-ACT). Free PSA is the other immunoreactive form present in serum. The third main form of PSA, bound to alpha-2-macroglobulin, cannot be detected by enzyme immunoassay tests.

PSA levels rise in prostatic pathologies such as benign prostatic hyperplasia (BPH) or prostate cancer (4). Testing for PSA and its evolution is useful for monitoring and controlling the efficacy of prostatic carcinoma therapy. The percentage of free PSA in serum is described as being significantly higher in patients with BPH than in patients with prostate cancer (5). Calculation of the percentage of free PSA, determined by dividing the free PSA (FPSA) concentration by that of total PSA (TPSA), has been suggested as a way of improving the differentiation of BPH and prostate cancer (2). The VIDAS FPSA assay is used to measure the free PSA concentration in order to calculate the free PSA percentage.

PRINCIPLE

The assay principle combines a two step enzyme immunoassay sandwich method with a final fluorescent detection (ELFA).

The Solid Phase Receptacle (SPR®), serves as the solid phase as well as the pipetting device for the assay. Reagents for the assay are ready-to-use and pre-dispensed in the sealed reagent strips.

All of the assay steps are performed automatically by the instrument. The sample is cycled in and out of the SPR several times. This operation enables the antibody fixed onto the interior wall of the SPR to capture the free fraction of the prostate specific antigen present in the sample. Unbound components are eliminated during the washing steps. The alkaline phosphatase-labeled antibody is then incubated in the SPR where it binds with the prostate specific antigen. Unbound conjugate is then eliminated during the washing steps.

During the final detection step, the substrate (4-Methyl-umbelliferyl phosphate) is cycled in and out of the SPR. The conjugate enzyme catalyzes the hydrolysis of this substrate into a fluorescent product (4-Methyl-umbelliferone) the fluorescence of which is measured at 450 nm. The intensity of the fluorescence is proportional to the concentration of prostate specific antigen free fraction present in the sample.

At the end of the assay, the results are automatically calculated by the instrument in relation to the calibration curve stored in memory, and then printed out.

CONTENT OF THE KIT (30 TESTS) - RECONSTITUTION OF REAGENTS:

30 FPSA strips	STR	Ready-to-use.
30 FPSA SPRs 1 x 30	SPR	Ready-to-use. Interior of SPRs coated with mouse monoclonal anti-free PSA specific immunoglobulins.
FPSA control 1 x 2 ml (lyophilized)	C1	Reconstitute with 2 ml of distilled water. Let stand for 30 minutes, then mix. Stable after reconstitution for 24 hours at 2-8°C or until expiration date on kit at -25 ± 6 °C. Five freeze/thaw cycles are possible. Human serum* + human PSA (free fraction) + preservatives. MLE data indicate the confidence interval in ng/mL ("Control C1 Dose Value Range").
FPSA calibrator 1 x 2.5 ml	S1	Ready-to-use. Bovine albumin + human PSA (free fraction) + 0.9 g/l sodium azide. MLE data indicate the concentration in ng/mL ("Calibrator (S1) Dose Value") and the confidence interval in "Relative Fluorescence Value" ("Calibrator (S1) RFV Range").

Specifications for the factory master data required to calibrate the test:

- MLE data (Master Lot Entry) provided in the kit, or
- MLE bar codes printed on the box label.
- 1 Package Insert provided in the kit or downloadable from www.biomerieux.com/techlib.
- * This product has been tested and shown to be negative for HBs antigen, antibodies to HIV1, HIV2 and HCV. However, since no existing test method can totally guarantee their absence, this product must be treated as potentially infectious. Therefore, usual safety procedures should be observed when handling.

The SPR

The interior of the SPR is coated during production with mouse monoclonal immunoglobulins directed against the free fraction of the prostate specific antigen. Each SPR is identified by the code "FPSA". Only remove the required number of SPRs from the pouch and carefully reseal the pouch after opening.

The strip

The strip consists of 10 wells covered with a labeled, foil seal. The label comprises a bar code which mainly indicates the assay code, kit lot number and expiration date. The foil of the first well is perforated to facilitate the introduction of the sample. The last well of each strip is a cuvette in which the fluorometric reading is performed. The wells in the center section of the strip contain the various reagents required for the assay.

Description of the FPSA strip:

Wells	Reagents				
1	Sample well.				
2 - 3 - 4	Empty wells.				
5	Conjugate: alkaline phosphatase-labeled monoclonal anti-PSA immunoglobulins (mouse) + Tris (0.1 mol/l, pH 6.5) + NaCl (0.1 mol/l) + calf serum (5%) + 0.9 g/l sodium azide (400 µl).				
6 - 7 - 9	Wash buffer: Tris (0.05 mol/l, pH 7.4) + NaCl (0.4 mol/l) + Tween (0.05%) + 1 g/l sodium azide (600 µl).				
8	Diluent: Tris (0.1 mol/l, pH 7.0) + NaCl (0.1 mol/l) + calf serum (5%) + sodium azide 0.9 g/l (400 μl).				
10	Cuvette with substrate: 4-Methyl-umbelliferyl phosphate (0.6 mmol/l) + diethanolamine (DEA*) (0.62 mol/l or 6.6%, pH 9.2) + 1 g/l sodium azide (300 µl).				

* Signal Word: DANGER



Hazard statement

H318: Causes serious eye damage.

Precautionary statement

P280 :Wear protective gloves/protective clothing/eye protection/face protection.

P305 + P351 + P338: IF IN EYES: Rinse cautiously with water for several minutes. Remove contact lenses, if present and easy to do. Continue rinsing.

For further information, refer to the Material Safety Data Sheet.

MATERIALS AND DISPOSABLES REQUIRED BUT NOT PROVIDED

- Pipette with disposable tip to dispense 2 ml and 200 µl.
- Powderless, disposable gloves.
- For other specific materials and disposables, please refer to the Instrument User's Manual.
- Instrument of the VIDAS family.

WARNINGS AND PRECAUTIONS

- For in vitro diagnostic use only.
- For professional use only.
- This kit contains products of human origin. No known analysis method can totally guarantee the absence of transmissible pathogenic agents. It is therefore recommended that these products be treated as potentially infectious and handled observing the usual safety precautions (see Laboratory biosafety manual - WHO - Geneva latest edition).
- This kit contains products of animal origin. Certified knowledge of the origin and/or sanitary state of the animals does not totally guarantee the absence of transmissible pathogenic agents. It is therefore recommended that these products be treated as potentially infectious and handled observing the usual safety precautions (do not ingest or inhale).

- Do not use the SPRs if the pouch is pierced.
- Do not use visibly deteriorated STRs (damaged foil or plastic).
- Do not use reagents after the expiration date indicated on the label.
- Do not mix reagents (or disposables) from different lots.
- Use **powderless** gloves, as powder has been reported to cause false results for certain enzyme immunoassay tests.
- Kit reagents contain sodium azide which can react with lead or copper plumbing to form explosive metal azides.
 If any liquid containing sodium azide is disposed of in the plumbing system, drains should be flushed with water to avoid build-up.
- The substrate in well 10 contains an irritant agent (6.6% diethanolamine). Refer to the hazard statements "H" and the precautionary statements "P" above.
- Spills should be wiped up thoroughly after treatment with liquid detergent or a solution of household bleach containing at least 0.5% sodium hypochlorite. See the User's Manual for cleaning spills on or in the instrument. Do not autoclave solutions containing bleach.
- The instrument should be regularly cleaned and decontaminated (see the User's Manual).

STORAGE CONDITIONS

Store the VIDAS FPSA kit at 2-8°C.

- Do not freeze SPRs, strips or calibrator.
- Store all unused reagents at 2-8°C.
- After opening the kit, check that the SPR pouch is correctly sealed and undamaged. If not, do not use the SPRs
- Carefully reseal the pouch with the desiccant inside after use to maintain stability of the SPRs and return the complete kit to 2-8°C.
- If stored according to the recommended conditions, all components are stable until the expiration date indicated on the label. Refer to the kit composition table for special storage conditions.

SPECIMENS

Specimen type and collection:

Serum or plasma (lithium heparin or EDTA).

As some collection tubes may contain substances which interfere with test results, it is recommended that each laboratory checks the compatibility of collection tubes used.

Samples containing impurities must be centrifuged before analysis.

None of the following factors have been found to significantly influence this assay.

- hemolysis (after spiking samples with hemoglobin: 0 to 300 μmol/l (monomer)),
- lipemia (after spiking samples with lipids: 0 to 10 mg/ml equivalent in triglycerides),
- bilirubinemia (after spiking samples with bilirubin: 0 to 470 μmol/l).

However, it is recommended not to use samples that are clearly hemolyzed, lipemic or icteric and, if possible, to collect a new sample.

Specimen stability:

From collection to assaying, samples (sera or plasma) can be stored at 2-8°C in stoppered tubes for up to 24 hours or at -25 \pm 6°C for up to 2 months. For better storage, avoid successive freezing and thawing. Some publications indicate that certain samples may lose their immunoreactivity during storage (6).

INSTRUCTIONS FOR USE

For complete instructions, see the User's Manual.

Reading Master lot data

Before each new lot of reagents is used, enter the specifications (or factory master data) into the instrument using the master lot entry (MLE) data.

If this operation is not performed **before initiating the tests**, the instrument will not be able to print results.

Note: the master lot data need only be entered once for each lot.

It is possible to enter MLE data **manually or automatically** depending on the instrument (refer to the User's Manual).

CALIBRATION

Calibration, using the calibrator provided in the kit, must be performed each time a new lot of reagents is opened, after the master lot data have been entered. Calibration should then be performed every 14 days. This operation provides instrument-specific calibration curves and compensates for possible minor variations in assay signal throughout the shelf-life of the kit.

The calibrator, identified by S1, must be tested in **duplicate** (see User's Manual). The calibrator value must be within the set RFV "Relative Fluorescence Value" range. If this is not the case, recalibrate.

VIDAS FPSA has been calibrated to the First International Standard for free Prostate-specific antigen (7). Depending on the dilution mode and the type of diluent used with the International Standard, a bias of up to 25% may be observed.

Procedure

- 1. Only remove the required reagents from the refrigerator and allow them to come to room temperature for at least 30 minutes.
- Use one "FPSA" strip and one "FPSA" SPR for each sample, control or calibrator to be tested. Make sure the storage pouch has been carefully resealed after the required SPRs have been removed.
- The test is identified by the "FPSA" code on the instrument. The calibrator must be identified by "S1", and tested in duplicate. If the control is to be tested, it should be identified by "C1".
- Mix the calibrator, control and samples using a vortextype mixer (for serum or plasma separated from the pellet).
- 5. For this test, the calibrator, control, and sample test portion is 200 μ l.
- Insert the "FPSA" SPRs and "FPSA" strips into the instrument. Check to make sure the color labels with the assay code on the SPRs and the Reagent Strips match.
- Initiate the assay as directed in the User's Manual. All the assay steps are performed automatically by the instrument.
- 8. Reclose the vials and return them to the required temperature after pipetting.
- The assay will be completed within approximately 60 minutes. After the assay is completed, remove the SPRs and strips from the instrument.
- 10. Dispose of the used SPRs and strips into an appropriate recipient.

RESULTS AND INTERPRETATION

Once the assay is completed, results are analyzed automatically by the computer. Fluorescence is measured twice in the Reagent Strip's reading cuvette for each sample tested. The first reading is a background reading of the substrate cuvette before the SPR is introduced into the substrate. The second reading is taken after incubating the substrate with the enzyme remaining on the interior of the SPR. The RFV (Relative Fluorescence Value) is calculated by subtracting the background reading from the final result. This calculation appears on the result sheet.

The results are automatically calculated by the instrument using calibration curves which are stored by the instrument (mathematical model: 4 parameter logistics model) the concentrations are expressed in ng/ml.

Samples with free PSA concentrations greater than 10 ng/ml should be retested after being diluted in the VIDAS TPSA kit diluent (ref. 30 428). If the dilution factor has not been entered when the Work List was created (see User's Manual), multiply the result by the dilution factor to obtain the sample concentration.

To determine the free PSA percentage, determine the Total PSA concentration for the same sample using VIDAS TPSA (ref. 30428) and the free PSA concentration using VIDAS FPSA (ref. 30440).

Calculate the percentage as follows:

FPSA concentration × 100

TFSA concentration

Interpretation of test results should be made taking into consideration the patient's history, and the results of any other tests performed.

QUALITY CONTROL

A control is included in each VIDAS FPSA kit. This control must be performed immediately after opening a new kit to ensure that reagent performance has not been altered. Each calibration must also be checked using this control. Results cannot be validated if the control value deviates from the expected values.

Note

It is the responsibility of the user to perform Quality Control in accordance with any local applicable regulations.

LIMITATIONS OF THE METHOD

Interference may be encountered with certain sera containing antibodies directed against reagent components. For this reason, assay results should be interpreted as part of a complete clinical profile.

The free PSA percentage can only be interpreted taking into consideration clinical data and information available from other diagnostic procedures.

Free PSA assays should not be performed on patients who have received a contrast medium within the previous 24 hours (8).

RANGE OF EXPECTED VALUES

A study was performed with the VIDAS FPSA and TPSA reagents using 227 samples (TPSA levels between 4 and 25 ng/ml) collected from patients with BPH (n=78) or with prostate cancer confirmed by biopsies (n=149). The statistical analysis of FPSA/TPSA ratios established at 18% the threshold associating the best sensitivity to the best specificity in the diagnosis of prostate cancer (FPSA/TPSA ratio < 18%).

		95% CI*
Sensitivity	73.15%	65.36 - 79.74%
Specificity	78.21%	67.61 - 86.05%
Positive predictive value	86.51%	79.29 - 91.48%
Negative predictive value	60.40%	50.45 - 69.55%

^{*95%} confidence interval

These figures are given as a guide; it is recommended that each laboratory establishes its own reference values from a rigorously selected population.

PERFORMANCE

Studies performed using VIDAS FPSA gave the following results:

Measurement range

The measurement range of the VIDAS FPSA kit extends up to 10 ng/ml.

Detection limit

Defined as the smallest concentration of PSA which is significantly different from the zero concentration with a probability of 95%: **0.05 ng/ml.**

Hook effect

No hook effect was found up to PSA free fraction concentrations of 50,000 ng/ml.

Precision

Within-run reproducibility:

Five samples were tested 30 times in the same run.

Sample	1	2	3	4	5
Mean (ng/ml)	0.63	1.59	3.10	5.66	9.69
CV %	3.92	3.29	2.17	3.33	3.15

Between-run reproducibility:

Five samples were tested in 29 different runs on the same VIDAS.

Sample	1	2	3	4	5
Mean (ng/ml)	0.67	1.69	3.45	6.15	9.79
CV %	5.20	4.42	4.04	3.85	3.84

Accuracy

Dilution test

Three samples were diluted in the VIDAS TPSA (ref. 30 428) kit diluent and tested singly in 3 runs. The ratio of the mean concentration measured over the expected concentration is expressed as a mean recovery percentage.

Sample no.	Dilution factor	Expected concentration (ng/ml)	Measured mean concentration (ng/ml)	Mean recovery percentage (%)
1	1/1	4.82	4.82	100.0
	1/2	2.41	2.42	100.6
	1/4	1.21	1.17	97.5
	1/8	0.60	0.59	98.4
	1/16	0.30	0.30	98.9
2	1/1	6.11	6.11	100.0
	1/2	3.05	3.17	103.8
	1/4	1.53	1.56	102.5
	1/8	0.76	0.81	106.3
	1/16	0.38	0.39	101.9
3	1/1	2.78	2.78	100.0
	1/2	1.39	1.33	95.6
	1/4	0.69	0.65	93.3
	1/8	0.35	0.34	98.2
	1/16	0.17	0.16	93.5

Comparison with other test methods

The concentration of PSA free fraction may vary in a sample determined using kits from different manufacturers, depending on the test methods used. If the test method is changed, laboratories should confirm the concentrations previously found. Correlation was established between the VIDAS FPSA kit (X) and a radioimmunoassay method (Y):

X = 1.15Y - 0.21 r = 0.990

CROSS REACTIVITY

(n = 71)

The immunological specificity of VIDAS FPSA was tested with a buffer spiked with PSA-ACT (bound PSA).

Spiking with	VIDAS FPSA	Recovery
PSA-ACT	assay	percentage
(ng/ml)	(ng/ml)	(%)
500	0.129	0.03

WASTE DISPOSAL

Dispose of used or unused reagents as well as any other contaminated disposable materials following procedures for infectious or potentially infectious products.

It is the responsibility of each laboratory to handle waste and effluents produced according to their nature and degree of hazardousness and to treat and dispose of them (or have them treated and disposed of) in accordance with any applicable regulations

LITERATURE REFERENCES

- OSTERLING JE. Prostate specific antigen: A critical assessment of the most useful tumor marker for adenocarcinoma of the prostate. J Urol 1991; <u>145</u>: 907-23.
- 2. CHRISTENSON A, LAURELL C-B, LILIJA H. Enzymatic activity of prostate-specific antigen and its reactions with extracellular serine proteinase inhibitor. Eur J Biochem 1990; 194: 755-63.
- ZHANG WM, LEINONEN J, KALKKINEN N, DOWELL B, STENMAN UH. Purification and characterization of different molecular forms of prostate-specific antigen in human seminal fluid. Clin Chem 1995; 41: 1567-1573.
- STAMEY TA. YANG N, HAY AR et al. Prostate specific antigen a serum marker for adenocarcinoma of the prostate. New Engl J Med 1987, 317: 909-16.
- STENMAN U-H, LEINONEN J, ZHANG W-H, Problems in the determination of prostate specific antigen. Eur J Clin Chem Clin Biochem 1996; 34: 735-40.
- PRICE C.P., et al. Pre- and post-analytical factors that may influence use of serum prostate specific antigen and its isoforms in a screening programme for prostate cancer. Ann Clin Biochem 2001; 38: 188-216.
- RAFFERTY B., RIGSBY P., ROSE. M., et al. Reference Reagents for Prostate-specific antigen (PSA): Establishment of the First International Standards for free PSA and PSA (90: 10). Clinical Chemistry, 2000, vol. 46, 1310-1317.
- WATANABE N. and al. In vitro effect of contrast during immunoradiometric assay for tumour-associated antigens. Nuclear Medicine Communication, 1998, 19, 63-70.

VIDAS[®] FPSA (FPSA) 09573 **H** - en - 2015/01

INDEX OF SYMBOLS

Symbol	Meaning	
REF	Catalog number	
IVD	In Vitro Diagnostic Medical Device	
	Manufacturer	
	Temperature limit	
	Use by date	
LOT	Batch code	
[]i	Consult Instructions for Use	
Σ	Contains sufficient for <n> tests</n>	
Date of manufacture		

WARRANTY

bioMérieux disclaims all warranties, express or implied, including any implied warranties of MERCHANTABILITY AND FITNESS FOR A PARTICULAR USE. bioMérieux shall not be liable for any incidental or consequential damages. IN NO EVENT SHALL BIOMERIEUX'S LIABLITY TO CUSTOMER UNDER ANY CLAIM EXCEED A REFUND OF THE AMOUNT PAID TO BIOMERIEUX FOR THE PRODUCT OR SERVICE WHICH IS THE SUBJECT OF THE CLAIM.

REVISION HISTORY

Change type categories:

N/A Not applicable (First publication)
Correction Correction of documentation anomalies

Technical change Addition, revision and/or removal of information related to the product Administrative Implementation of non-technical changes noticeable to the user

Note: Minor typographical, grammar, and formatting changes are not included in the

revision history.

	Release date	Part Number	Change Type	Change Summary
2015/01		Administrative	INDEX OF SYMBOLS REVISION HISTORY	
	2015/01	09573H	Technical	CONTENT OF THE KIT (30 TESTS) – RECONSTITUTION OF REAGENTS WARNINGS AND PRECAUTIONS INSTRUCTIONS FOR USE

BIOMERIEUX, the blue logo, SPR and VIDAS are used, pending and/or registered trademarks belonging to bioMérieux, or one of its subsidiaries, or one of its companies.

Any other name or trademark is the property of its respective owner.





